



# Ultrasound Tables

• Ultrasound Tables • Echocardiography Tables • Vascular Tables



**NEW** Adjustable Headrest Extension

**NEW** Motor-controlled ergonomic Fowler back with infinite adjustments through 80°

**NEW** Durable Hand Controller

Paper dispenser

**NEW** Lever to return cardiac cushion to level

**NEW** Industry's largest cardiac scanning cutout

Sonographer cutout

**NEW** Articulating Scanning Arm Board<sup>1</sup>

**NEW** Lever to drop down cardiac cushion

**NEW** Flush-mounted Side Rails fold beneath table for unencumbered access<sup>1</sup>

**NEW** Retractable Stirrups with contoured heel cup for ergonomic support<sup>1</sup>

All motorized tables are ETL product listed to the UL 60601-1 Medical Electrical Equipment Safety Standard.

<sup>1</sup>Standard on Sound Pro, Optional on all other tables

Foot Controller (Standard or optional on select tables, not shown)

Central Floor-Locking System (Standard on select tables, not shown)

[www.biodex.com/ultrasound](http://www.biodex.com/ultrasound)

# BIODEX **Ultrasound**Tables

Say goodbye to yesterday,  
and hello to tomorrow...

How do you improve the number one selling Ultrasound Tables in the world?

By simply making them better... make the scanning  
experience more comfortable for the patient, and more  
ergonomically correct for the sonographer.

And that is just what Biodex has done.

With customer satisfaction driving the decision, Biodex has applied  
hands-on suggestion to the redesign.

The result... feature-rich tables that meet the concerns  
of today's sonographer, at a cost you can afford.

Extra-wide tabletop to  
accommodate bariatric population

Fold down leg section

**WARNING**  
2000 LBS. LOAD CAPACITY  
Prohibited Use:  
Do not use for other body conditions.  
Restrictions:  
Do not place or perform any other body condition  
exercises which are prohibited by this table.



# SoundPro™

## Suitable Procedures:

- Abdomen
- Renal
- Aorta
- Pelvis
- Endovaginal
- OB
- Hysterosonography
- Thyroid
- Testicular
- Breast
- Musculoskeletal
- Carotid
- DVT/Lower Extremity Venous
- Lower Extremity Arterial
- Venous Reflux
- Upper Extremity Venous
- Upper Extremity Arterial
- Renal Artery
- Renal Transplant
- Abdominal Doppler
- SMA
- Echocardiography
- Stress Echo
- Transesophageal Echo

*Accommodates*  
**BARIATRIC**  
*Patients*



## -Enhanced ergonomics - accommodates ultrasound and echocardiography procedures.

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra wide top with 500-lb patient capacity, suitable for bariatric patients
- Cardiac drop-down cushion, allows release and return from either side of table
- Industry's largest cardiac scanning cutout – 9.5" w x 15.75" l
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment
- Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller activates table positions
- Allows seated or standing work positions
- Redesigned side rails mount flush to table for unencumbered patient access
- Enhanced stirrups, designed to cup at the heel for added comfort
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- Foot Controller is ideal for when hands are busy with other functions
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.

### Optional

- Headrest extension drops down 15 degrees to allow for carotid and thyroid studies

*All table options and accessories are featured on pages 18-19*

[www.biodex.com/soundpro](http://www.biodex.com/soundpro)

The Sound Pro™ Combination Table brings together the imaging features of both the Ultra Pro and Echo Pro tables. Fully equipped for OB/GYN, general ultrasound and echocardiography procedures, the extra wide tabletop is designed with advanced ergonomics to provide



maximum comfort for both patient and sonographer. Designed to conform to natural body extension and sitting position, the motorized Fowler back is infinitely adjustable to 80 degrees via hand or foot control.

The table is so

accommodating, it reduces the total amount of time required to achieve a quality image. Stress patients are off the treadmill and quickly onto the table in imaging position.



The new cardiac drop-down cushion releases from either side of the table, providing open access to the patient's left thorax area for an unencumbered apical approach. Flush-mounted side rails enable closer access to the

patient, without sonographer contortion. When not in use, rails fold beneath the table, eliminating obstruction.

Patient comfort is assured even during lengthy procedures with our new softer feel, Naugahyde® antimicrobial mattress cover, which also provides protection against bacteria. Standard accessories available to expand the function of your table include ergonomic stirrups that cup at the heel, an articulating scanning arm board (adjustable from 0 to 130 degrees) and a foot controller for when hands are busy with other functions. The NEW, optional Headrest Extension is removable and adjustable, providing convenient access and patient comfort during carotid and thyroid procedures.



*Biodex continually uses customer input as we redesign our tables to enhance ergonomics for the comfort and safety of the sonographer.*

### Sound Pro™ Combination Table

**058-710** Table, Combination, Sound Pro™, 115 VAC

*Includes: hand and foot controllers, folding side rails, articulating scanning arm board, retractable stirrups, IV pole and paper dispenser.*

**058-715** Table, Combination, Sound Pro™, 230 VAC

*Includes: hand and foot controllers, folding side rails, articulating scanning arm board, retractable stirrups, IV pole and paper dispenser.*

### Related:

**058-704** Chair, Ergonomic, Sonography

### Optional:

**058-738** Extension, Headrest

*Specifications and models listed on pages 20-21*

**BIODEX**



# UltraPro™

## Suitable Procedures:

- Abdomen
- Renal
- Aorta
- Pelvis
- Endovaginal\*\*
- OB
- Hysterosonography
- Thyroid
- Testicular
- Breast
- Musculoskeletal
- Carotid
- DVT/Lower Extremity Venous
- Lower Extremity Arterial
- Upper Extremity Venous\*
- Renal Artery
- Renal Transplant
- Abdominal Doppler
- SMA
- Transesophageal Echo

\* Requires articulating scanning arm board.

\*\*Requires optional stirrups.

*Accommodates*  
**BARIATRIC**  
*Patients*



## Comfort focused design for quality imaging

- Ergonomic design for enhanced patient and sonographer comfort
- 500-lb patient capacity, suitable for bariatric patients
- Drop-down, contoured leg section
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment
- Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller activates table positions
- Allows seated or standing work positions

### Optional

- Foot Controller is ideal for when hands are busy with other functions
- Redesigned side rails mount flush to table for unencumbered patient access
- Enhanced stirrups, designed to cup at the heel for added comfort
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- Headrest extension for easier thyroid and carotid access
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.

*All table options and accessories are featured on pages 18-19*

The Ultra Pro™ Ultrasound Table is the professional's choice when performing general ultrasound procedures, including OB/GYN.

Now there is even more reason to choose an Ultra Pro Ultrasound Table. Enhancements to the ergonomics reduce sonographer fatigue and require less patient repositioning, making it easy to achieve quality images in less time.



The redesigned Fowler back which provides the utmost patient comfort, an enhanced central floor-locking system offering easier access from around the table and an antimicrobial vinyl mattress ensuring safety against infection are all standard on the Ultra Pro.



Customization is made possible with optional items including flush mounted side rails allowing complete patient access and safety, retractable stirrups with enhanced contouring and extra-long extension for taller patients, an articulating scanning arm board and a headrest extension that drops down 15 degrees to allow for carotid studies.



*The motorized Fowler back has minimal patient shift of no more than 2-3 inches.*

### Ultra Pro™ Combination Table

**058-720** Table, Ultrasound, Ultra Pro™, 115 VAC  
*Includes hand controller.*

**058-725** Table, Ultrasound, Ultra Pro™, 230 VAC  
*Includes hand controller.*

### Related:

**058-704** Chair, Ergonomic, Sonography

### Optional:

- 058-740** Foot Controller
- 058-633** Side Rails, Folding
- 058-652** Stirrups, Retractable
- 058-736** Arm Board, Articulating Scanning
- 058-737** IV Pole
- 058-611** Dispenser, Paper
- 058-738** Extension, Headrest

*Specifications and models listed on pages 20-21*



# EconoUltrasound

## Suitable Procedures:

- Abdomen
- Renal
- Aorta
- Pelvis
- Endovaginal\*\*
- OB
- Hysterosonography
- Thyroid
- Testicular
- Breast
- Musculoskeletal
- Carotid
- DVT/Lower Extremity Venous
- Lower Extremity Arterial
- Venous Reflux
- Upper Extremity Venous
- Upper Extremity Arterial\*
- Renal Artery
- Renal Transplant
- Abdominal Doppler
- SMA
- Transesophageal Echo

\* Requires articulating scanning arm board.

\*\*Requires optional stirrups.

*Accommodates*  
**BARIATRIC**  
*Patients*





# Ultrasound Table

## Economically priced, rich with features.

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- 500-lb patient capacity, suitable for bariatric patients
- Drop-down, contoured leg section
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height and Fowler positioning
- Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Durable hand controller activates table positions
- Allows seated or standing work positions

### Optional

- Redesigned side rails mount flush to table for unencumbered patient access
- Enhanced stirrups, designed to cup at the heel for added comfort
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- Headrest extension for easier thyroid and carotid access
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.

*All table options and accessories are featured on pages 18-19*

The Econo Ultrasound Table is designed on the same platform as the Ultra Pro with similar features and functions. To make it a more cost-effective alternative it is offered without Trendelenburg positioning or central floor-locking system.

The same enhanced ergonomics have been incorporated into the design of the Econo Ultrasound Table, making it easy for a sonographer to achieve quality images without discomfort to the patient or sonographer.



Motorized height and Fowler positioning which accommodates natural body extension and sitting position, a height range as low as 23" to allow wheelchair transfers and a new hand controller are just some of the features that make the imaging experience easy with the Econo Ultrasound Table.



The Econo Table can be customized to suit your department needs with an assortment of table options. Flush-mounted side rails enable closer patient access without sonographer contortion. A vascular scanning arm board, adjustable from 0-130 degrees and a headrest extension that drops down 15 degrees to allow for carotid studies are just some of the items which can be added to the table.



*The ergonomics of  
an ultrasound table  
are crucial to lowering  
the risk of incurring  
musculoskeletal injuries.*

### Econo Ultrasound Table

**058-726** Table, Ultrasound, Econo, 115 VAC  
*Includes hand controller.*

**058-727** Table, Ultrasound, Econo, 230 VAC  
*Includes hand controller.*

### Related:

**058-704** Chair, Ergonomic, Sonography

### Optional:

**058-633** Side Rails, Folding  
**058-652** Stirrups, Retractable  
**058-736** Arm Board, Articulating Scanning  
**058-737** IV Pole  
**058-611** Dispenser, Paper  
**058-738** Extension, Headrest

*Specifications and models listed on pages 20-21*

# EchoPro™

## Suitable Procedures:

- Echocardiography
- Stress Echo
- Transesophageal Echo
- Abdomen
- Renal
- Aorta
- Pelvis
- Thyroid
- Testicular
- Breast
- Musculoskeletal
- Carotid
- Upper Extremity Venous\*
- Lower Extremity Arterial
- Upper Extremity Arterial\*
- Renal Artery
- Renal Transplant
- Abdominal Doppler
- SMA

*Accommodates*  
**BARIATRIC**  
*Patients*



*\* Requires articulating scanning arm board.*



# Echocardiography Table

## Improved sonographer comfort - closer patient access

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra wide top with 500-lb patient capacity, suitable for bariatric patients
- Cardiac drop-down cushion, allows release and return from either side of table
- Industry's largest cardiac scanning cutout – 9.5" w x 15.75" l
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment
- Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller activates table positions
- Allows seated or standing work positions


### Optional

- Foot Controller is ideal for when hands are busy with other functions
- Redesigned side rails mount flush to table for unencumbered patient access
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
- Headrest extension drops down 15 degrees to allow for carotid and thyroid studies

*All table options and accessories are featured on pages 18-19*

The Echo Pro™ Echocardiography Table has been re-engineered with the sonographer in mind. The new ergonomic design provides a barrier-free imaging experience. The result is less sonographer injury and better quality images.

Redesigned features make it easier for the sonographer to access the patient from all sides. The new cardiac scanning cushion now drops down and releases from either side of the table with nothing protruding. The sonographer has open access to the left thorax area for an unobstructed apical approach. New side rails are flush to the table and fold completely beneath the table out of the sonographer's way, eliminating potential obstruction.



Patient comfort is assured even during lengthy procedures with our new softer feel, Naugahyde® antimicrobial mattress cover, which also provides protection against bacteria.

Add options such as an Articulating Scanning Arm Board or



*The cardiac scanning cushion has been redesigned to drop down and release from either side of the table, with nothing protruding.*

### Echo Pro™ Echocardiography Table

**058-700** Table, Echocardiography, Echo Pro™, 115 VAC  
*Includes hand controller.*

**058-705** Table, Echocardiography, Echo Pro™, 230 VAC  
*Includes hand controller.*

### Related:

**058-704** Chair, Ergonomic, Sonography

### Optional:

**058-740** Foot Controller  
**058-633** Side Rails, Folding  
**058-736** Arm Board, Articulating Scanning  
**058-737** IV Pole  
**058-611** Dispenser, Paper  
**058-738** Extension, Headrest

*Specifications and models listed on pages 20-21*

# EconoEchocardiography

*Accommodates*  
**BARIATRIC**  
*Patients*

## Suitable Procedures:

- Echocardiography
- Stress Echo
- Transesophageal Echo
- Abdomen
- Renal
- Aorta
- Pelvis
- Thyroid
- Testicular
- Breast
- Musculoskeletal
- Carotid
- Upper Extremity Venous\*
- Upper Extremity Arterial\*
- Renal Artery
- Renal Transplant
- Abdominal Doppler
- SMA

\* Requires articulating scanning arm board.





# Echocardiography Table

## A feature rich table, economically priced

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra wide top with 500-lb patient capacity, suitable for bariatric patients
- Cardiac drop-down cushion, allows release and return from either side of table
- Industry's largest cardiac scanning cutout – 9.5" w x 15.75" l
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height and Fowler positioning
- Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Durable hand controller activates table positions
- Allows seated or standing work positions

### Optional

- Redesigned side rails mount flush to table for unencumbered patient access
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
- Headrest extension drops down 15 degrees to allow for carotid and thyroid studies

*All table options and accessories are featured on pages 18-19*

The Econo Echocardiography table is based on the same platform as the Echo Pro. It has the same ergonomic features, keeping the safety of the sonographer in mind. To make it a more cost-effective alternative, it is offered without Trendelenburg positioning or the central floor-locking system.



The barrier between sonographer and patient is eliminated with the Econo Echocardiography Table, featuring a new ergonomic design that helps prevent musculoskeletal injury.

The new design includes a cardiac drop-down cushion that provides the largest scanning access of any table, with the ability to release and return from either side of the table. Open access is provided to the patient's left thorax area for an unobstructed apical approach.



The Econo Echocardiography Table has many optional items to further improve sonographer comfort and safety. Choose what best suit your department's specific needs.

Redesigned side rails for closer patient access, a vascular scanning arm board, adjustable from 0 to 130 degrees and a headrest extension for easier thyroid and carotid access are all available.



*Biodex Ultrasound table enhancements - inspired by customer feedback.*

### Econo Echocardiography Table

**058-701** Table, Echocardiography, Econo, 115 VAC  
*Includes hand controller.*

**058-706** Table, Echocardiography, Econo, 230 VAC  
*Includes hand controller.*

### Related:

**058-704** Chair, Ergonomic, Sonography

### Optional:

- 058-633** Side Rails, Folding
- 058-736** Arm Board, Articulating Scanning
- 058-737** IV Pole
- 058-611** Dispenser, Paper
- 058-738** Extension, Headrest

*Specifications and models listed on pages 20-21*

# Echo/VascPro™

*Accommodates*  
**BARIATRIC**  
*Patients*

## Suitable Procedures:

- Echocardiography
- Stress Echo
- Transesophageal Echo

## Lower Extremity

- Venous Ablation
- Venous Insufficiency
- Venous Reflux
- Venous Mapping and Doppler (DVT)
- Arterial Scan/Doppler (PAD)

## Upper Extremity

- Arterial\*
- Testicular
- Musculoskeletal
- Carotid
- Carotid Duplex
- Transcranial
- Abdomen
- Aorta
- Thyroid
- Breast
- Renal
- Renal Transplant
- Renal Artery
- Venous Access Placement and Assessment\*
- SMA



Features 30° Reverse Trendelenburg

\* Requires articulating scanning arm board.



# Echocardiography and Vascular Table

## Achieve quality vascular and echocardiography images

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra wide top with 500-lb patient capacity, suitable for bariatric patients
- Cardiac drop-down cushion, allows release and return from either side of table
- Industry's largest cardiac scanning cutout – 9.5" w x 15.75" l
- Reverse Trendelenburg to 30 degrees
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment
- Improved foot support with comfortable access to calves without crouching
- Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller activates table positions
- Allows seated or standing work positions

### Optional

- Foot Controller is ideal for when hands are busy with other functions
- Redesigned side rails mount flush to table for unencumbered patient access
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
- Headrest extension drops down 15 degrees to allow for carotid and thyroid studies

All table options and accessories are featured on pages 18-19

[www.biodex.com/echovascpro](http://www.biodex.com/echovascpro)

The Echo/Vasc Pro™ Vascular Echocardiography Table is the logical choice for performing vascular scanning and



echocardiography. It has been redesigned with many features that help obtain the best quality image and keep the sonographer safe from injury.

The new Echo/Vasc Pro has

30 degree reverse Trendelenburg positioning, expanding lower extremity veins for maximum vessel dilation. Enhanced foot support provides the sonographer comfortable

access to a patient's calves without crouching. Body straps stabilize the patient and ensure safety throughout the entire procedure. The new design includes a cardiac drop-down cushion, with the ability to release and return from either side of the table, providing ample scanning access for an unobstructed apical approach.



As with all Biodex ultrasound tables the Echo/Vasc Pro has enhanced Fowler positioning to accommodate natural body extension and sitting position without slide and shift, redesigned optional side rails for easy

access to the patient, an optional Articulating Scanning Arm Board adjustable from 0 to 130 degrees and a new soft antimicrobial vinyl cover.



*One of the most important elements in ultrasound table ergonomic design is to optimize the sonographer's reach. It's essential for the sonographer to get as close to the patient as possible when performing a scanning procedure.*

### Echo/Vasc Pro™ Table

**058-702** Table, Vascular Echocardiography, Echo/Vasc Pro™, 115 VAC  
*Includes hand controller and foot support.*

**058-707** Table, Vascular Echocardiography, Echo/Vasc Pro™, 230 VAC  
*Includes hand controller and foot support.*

### Related:

**058-704** Chair, Ergonomic, Sonography

### Optional:

**058-740** Foot Controller  
**058-633** Side Rails, Folding  
**058-736** Arm Board, Articulating Scanning  
**058-737** IV Pole  
**058-611** Dispenser, Paper  
**058-738** Extension, Headrest

*Specifications and models listed on pages 20-21*

# BIODEX

# VascPro™

## Suitable Procedures:

### Lower Extremity

- Venous Insufficiency
- Venous Ablation
- Venous Reflux
- Venous Mapping and Doppler (DVT)
- Arterial Scan/Doppler (PAD)

### Upper Extremity

- Arterial\*
- Testicular
- Musculoskeletal
- Carotid
- Carotid Duplex
- Transcranial
- Abdomen
- Aorta
- Thyroid
- Breast
- Renal
- Renal Transplant
- Renal Artery
- Venous Access Placement and Assessment\*
- SMA

\* Requires articulating scanning arm board.

*Accommodates*  
**BARIATRIC**  
*Patients*



Features 30° Reverse Trendelenburg



# Vascular Ultrasound Table

## The ergonomic choice for long-leg vascular studies.

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra wide top with 500-lb patient capacity, suitable for bariatric patients
- Reverse Trendelenburg to 30 degrees
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment
- Improved foot support with comfortable access to calves without crouching
- Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller
- Allows seated or standing work positions

### Optional

- Foot Controller is ideal for when hands are busy with other functions
- Redesigned side rails mount flush to table for unencumbered patient access
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
- Headrest extension drops down 15 degrees to allow for carotid and thyroid studies

*All table options and accessories are featured on pages 18-19.*

The new Vasc Pro™ Vascular Ultrasound Table is designed with ergonomics in mind. It addresses sonographer comfort



and safety and proper patient positioning, virtually eliminating repetitive-strain injuries from vascular scanning.

Sonographer comfort and safety are accomplished through new features such as

the enhanced foot support allowing easier access to the patient calves without crouching and optional side rails that enable closer access to the patient without sonographer contortion.

Fowler positioning to 80 degrees, which accommodates natural body extension and sitting position without slide or shift, an optional Articulating Scanning Arm Board adjustable from 0 to 130-degrees locking every 10 degrees, the thick antimicrobial mattress and table adjustability all create a comfortable imaging experience for the patient.



*Optional side rails tuck underneath the table when not in use to allow clear access to the tabletop.*

### Vasc Pro™ Vascular Ultrasound Table

**058-732** Table, Vascular Ultrasound, Vasc Pro™, 115 VAC  
*Includes hand controller and foot support.*

**058-733** Table, Vascular Ultrasound, Vasc Pro™, 230 VAC  
*Includes hand controller and foot support.*

### Related:

**058-704** Chair, Ergonomic, Sonography

### Optional:

**058-740** Foot Controller  
**058-633** Side Rails, Folding  
**058-736** Arm Board, Articulating Scanning  
**058-737** IV Pole  
**058-611** Dispenser, Paper  
**058-738** Extension, Headrest

*Specifications and models listed on pages 20-21*

# Table Accessories

## ERGONOMIC Sonography Chair

*Modern, ergonomic design - Sonography Chair features a foot ring for additional support and positioning.*



Modern design and manufactured with recycled materials, the Sonography Chair presents an ergonomic enhancement when scanning patients. Multiple adjustments for seat height, back height and backwards tilt tension adjustment with integrated seat slider allow for a variety of comfortable postures. A foot ring is standard for additional support and positioning. Color is graphite, matching the 700 series of Biodex Ultrasound Tables.

**058-704** Chair, Ergonomic, Sonography  
*Includes Foot Ring*

Visit OSHA's website at [www.osha.gov](http://www.osha.gov) on "Sonography: Positioning Patients and Equipment" which discusses not only the sonography chair but also exam tables. Biodex meets all of the sonographer's needs!

## Foot Controller

*Hands busy? Use Foot Control.*



The ergonomically designed Foot Controller offers the freedom to activate height, Trendelenburg motions and Fowler positioning from locations away from the table, or when hands are busy with other functions. The controller is manufactured from a metal casting for extreme durability.

**058-740** Foot Controller  
*For 058-700, 058-702, 058-720 and 058-732 Ultrasound Tables.*

## IV Pole

*Easily attaches for quick convenience*



**058-737** IV Pole

The optional IV Pole mounts in a bracket at the head end of the table. Constructed of stainless steel, this IV Pole is height adjustable.

## Articulating Scanning Arm Board

*Adjustable from 0° - 130°, locking every 10°*



For the gamut of upper extremity vascular procedures, the optional arm board is used to support the patient's arm and can be mounted to either side of the table. The arm board is adjustable from 0 to 130 degrees, locking every 10 degrees. It is required to perform upper extremity arterial and upper extremity venous procedures.

**058-736** Arm Board, Articulating Scanning



# Expand the function of your table.

## Paper Dispenser

*Includes Cutter Strap*



The optional Paper Dispenser mounts to the head end of the table. The dispenser is manufactured to hold a roll of hygienic table paper. Includes cutter strap.

**058-611** Dispenser, Paper  
*with Cutter Strap*

**056-612** Table Paper, Hygienic, 12 rolls  
*18" w x 225' l (46 cm x 69 m)*

## Headrest Extension

*Ideal for thyroid and carotid studies*



The NEW optional Headrest Extension is removable and adjustable, providing convenient access and patient comfort. Extension drops down 15° for easy access to perform thyroid and carotid studies. It's versatility allows you to extend the length of the table to accommodate taller patients.

**058-738** Extension, Headrest

## Folding Side Rails

*Safety and Convenience*



Fold-away side rails provide added safety without obstructing patient access. When stowed, the side rails fold out of the way, allowing the sonographer total patient access and comfortable scanning positioning.

**058-633** Folding Side Rails

## Retractable Stirrups

*Contoured heel cup for comfortable support*



Retractable Stirrups expand function of the table to include endovaginal, hysterosonography, testicular, OB and more. Bottom cushion of the table drops out of the way for unobstructed access.

**058-652** Stirrups, Retractable

## The Table Adjusts...

### ...not you, not the patient.

Biodex Ultrasound Tables are "Sonographer Friendly". They allow the sonographer to adjust the table to ergonomically correct positions, thus avoiding sonographer pain or injury. Biodex tables can shorten exam time and improve results. Whether you are looking for an ultrasound table for OB/GYN, echocardiography or vascular applications, Biodex is the answer.

Biodex Medical Imaging Tables are engineered and manufactured to only one level of quality... the best. Dependability and exceptional performance allow you to forget about the product and focus on the patient - that's your "Clinical Advantage." Combine that with competitive prices and Biodex becomes the obvious choice.

*Options are installed by the customer.*

# BIODEX

# Which Table is right for you?

## Sound Pro Combination Table



### SPECIFICATIONS:

#### Dimensions:

Overall: 71" l x 30" w (180.3 x 76.2 cm); 35" w (88.9 cm) with side rails  
 Cardiac Scanning Cutout: 9.5 w x 15.75 l (24.1 x 40 cm)  
 Drop-down cushion releases from either side of the table.  
 Sonographer Cutout: 8" w x 14.75" l (20.3 x 37.4 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.  
 Drop-Down Leg Section: 12.6" l (32 cm)

#### Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)  
 Trendelenburg: 0° to ±15°  
 Fowler Back: 0° to 80° infinitely adjustable

#### Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions  
 Foot Control: Activates height, Trendelenburg motions and Fowler positioning  
 Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

**Tabletop:** Three primary sections: torso section features sonographer and patient cutouts with Fowler positioning; center section remains fixed and leg section drops down to 40° and 80° for stirrup access  
**Drop-Down Cardiac Scanning Cushion:** Release from either side of the table  
**Adjustable Sonographer Cushion:** drops down or folds up for maximum access to the patient

**Mattress:** Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria  
**Upholstery Color:** Graphite  
**Patient Restraints:** One body strap  
**Finish:** Powder coat  
**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.  
**Weight:** 350 lb (159 kg)  
**Power:** 115 VAC or 230 VAC  
**Warranty:** **Two-year**  
**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN 60601-1 standards.



**058-710** Table, Combination, Sound Pro, 115 VAC  
*Includes: band and foot controllers, folding side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.*

**058-715** Table, Combination, Sound Pro, 230 VAC  
*Includes: band and foot controllers, folding side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.*

## Ultra Pro Ultrasound Table



### SPECIFICATIONS:

#### Dimensions:

Overall: 70" l x 30" w (177.8 x 76.2 cm); 35" w (88.9 cm) with optional side rails  
 Drop-Down Leg Section: 12.6" l (32 cm)

#### Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)  
 Trendelenburg: 0° to ±15°  
 Fowler Back: 0° to 80° infinitely adjustable

#### Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions  
 Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning  
 Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

**Tabletop:** Three primary sections: torso section with Fowler positioning, center section remains fixed and leg section drops down to 40° and 80° for stirrup access

**Mattress:** Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria  
**Upholstery Color:** Graphite  
**Patient Restraints:** Two body straps  
**Finish:** Powder coat  
**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.  
**Weight:** 335 lb (152 kg)  
**Power:** 115 VAC or 230 VAC  
**Warranty:** **Two-year**  
**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN 60601-1 standards.



**058-720** Table, Ultrasound, Ultra Pro, 115 VAC  
*Includes Hand Controller.*

**058-725** Table, Ultrasound, Ultra Pro, 230 VAC  
*Includes Hand Controller.*

## Econo Ultrasound Table



### SPECIFICATIONS:

#### Dimensions:

Overall: 70" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails  
 Drop-Down Leg Section: 12.6" l (32 cm)

#### Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)  
 Fowler Back: 0° to 80° infinitely adjustable

#### Controls:

Hand Control: Activates height and Fowler back  
 Wheels: 5" (12.7 cm) Individual locking swivel casters

**Tabletop:** Three primary sections: torso section with Fowler positioning, center section remains fixed and leg section drops down to 40° and 80° for stirrup access

**Mattress:** Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria  
**Upholstery Color:** Graphite  
**Patient Restraints:** Two body straps  
**Finish:** Powder coat  
**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.  
**Weight:** 335 lb (152 kg)  
**Power:** 115 VAC or 230 VAC  
**Warranty:** **Two years**  
**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN 60601-1 standards.



**058-726** Table, Ultrasound, Econo, 115 VAC  
*Includes band controller.*

**058-727** Table, Ultrasound, Econo, 230 VAC  
*Includes band controller.*

See pages 18-19 for Table Accessories, available separately.



## Echo Pro Echocardiography Table



### SPECIFICATIONS:

#### Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm)  
with optional side rails  
Cardiac Scanning Cutout: 9.5 w x 15.75 l (24.1 x 40 cm)  
Drop-down cushion releases from either side of the table.  
Sonographer Cutout: 8" w x 14.75" l (20.3 x 37.4 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.

#### Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)  
Trendelenburg: 0° to ±15°  
Fowler Back: 0° to 80° infinitely adjustable

#### Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions  
Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

**Tabletop:** Two primary sections: torso section: features sonographer and patient cutouts with Fowler positioning, and lower body section remains fixed

**Mattress:** Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria  
**Upholstery Color:** Graphite  
**Patient Restraints:** Two body straps  
**Finish:** Powder coat  
**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.  
**Weight:** 335 lb (152 kg)  
**Power:** 115 VAC or 230 VAC  
**Warranty:** Two-year  
**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN60601-1 standards.



**058-700** Echocardiography, Echo Pro™, 115 VAC  
*Includes hand controller.*

**058-705** Table, Echocardiography, Echo Pro™, 230 VAC  
*Includes hand controller.*

## Econo Echocardiography Table



### SPECIFICATIONS:

#### Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm)  
with optional side rails  
Cardiac Scanning Cutout: 9.5 w x 15.75 l (24.1 x 40 cm)  
Drop-down cushion releases from either side of the table.  
Sonographer Cutout: 8" w x 14.75" l (20.3 x 37.4 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.

#### Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)  
Fowler Back: 0° to 80° infinitely adjustable

#### Controls:

Hand Control: Activates height and Fowler back  
Wheels: 5" (12.7 cm) Individual locking swivel casters

**Tabletop:** Two primary sections: torso section: features sonographer and patient cutouts with Fowler positioning, and lower body section remains fixed

**Mattress:** Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria  
**Upholstery Color:** Graphite  
**Patient Restraints:** Two body straps  
**Finish:** Powder coat  
**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.  
**Weight:** 335 lb (152 kg)  
**Power:** 115 VAC or 230 VAC  
**Warranty:** Two-year  
**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN60601-1 standards.



**058-701** Table, Echocardiography, Econo, 115 VAC  
*Includes hand controller.*

**058-706** Table, Echocardiography, Econo, 230 VAC  
*Includes hand controller.*

## Echo/Vasc Pro Vascular Ultrasound Table



### SPECIFICATIONS:

#### Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm)  
with optional side rails  
Cardiac Scanning Cutout: 9.5 w x 15.75 l (24.1 x 40 cm)  
Drop-down cushion releases from either side of the table.  
Sonographer Cutout: 8" w x 14.75" l (20.3 x 37.4 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.

#### Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)  
Trendelenburg: 0° to 15°  
Reverse Trendelenburg: 0° to 30°  
Fowler Back: 0° to 80° infinitely adjustable

#### Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions  
Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

**Tabletop:** Two primary sections: torso section features sonographer and patient cutouts with Fowler positioning, and lower body section remains fixed

**Mattress:** Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria  
**Upholstery Color:** Graphite  
**Patient Restraints:** Two body straps  
**Finish:** Powder coat  
**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.  
**Weight:** 335 lb (152 kg)  
**Power:** 115 VAC or 230 VAC  
**Warranty:** Two-year  
**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN60601-1 standards.



**058-702** Table, Vascular Ultrasound, Echo/Vasc Pro, 115 VAC  
*Includes hand controller and foot support.*

**058-707** Table, Vascular Ultrasound, Echo/Vasc Pro, 230 VAC  
*Includes hand controller and foot support.*

## Vasc Pro Vascular Ultrasound Table



### SPECIFICATIONS:

#### Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm)  
with optional side rails

#### Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)  
Trendelenburg: 0° to 15°  
Reverse Trendelenburg: 0° to 30°  
Fowler Back: 0° to 80° infinitely adjustable

#### Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions  
Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

**Tabletop:** Two primary sections: torso section with Fowler positioning; lower body section remains fixed

**Mattress:** Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria  
**Upholstery Color:** Graphite  
**Patient Restraints:** Three body straps  
**Finish:** Powder coat  
**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.  
**Weight:** 335 lb (152 kg)  
**Power:** 115 VAC or 230 VAC  
**Warranty:** Two-year  
**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN60601-1 standards.



**058-732** Table, Vascular Ultrasound, Vasc Pro™, 115 VAC  
*Includes hand controller and foot support.*

**058-733** Table, Vascular Ultrasound, Vasc Pro™, 230 VAC  
*Includes hand controller and foot support.*

## ■ PEER PERSPECTIVE

*“We find that this table is easy to use and we are pleased with its performance.”*

### Sound Pro Table



“Some of the numerous features that we find exceptional on the Sound Pro, are the pivoting side rails, the retractable stirrups, and the cardiac cutouts. The construction of the side rails allows the sonographer the ability to be closer to the patient without any obstruction of the side rails hitting our knees.

Furthermore, the cardiac cutout allows the sonographer the ease to scan apical views of the heart. Lastly, we like the fact that the stirrups are retractable and hidden in the table when we do not require them. We find that this table is easy to use and we are pleased with its performance. We look forward to the next table we purchase.”

*- Heather Hibner, RT, RDMS, RDCS,  
Lead Sonographer, Madera Community Hospital, Madera, CA*

*“Because the table is so accommodating, it reduces the time needed to get a quality image.”*

### Echo Pro Echocardiography Table

“Our sonographers really like the Biodex Echo Pro Table. They love the electronic height adjustment. It allows the table to be lowered to a level where even patients with limited mobility can get on it easily. Once the patient is on the table, you can then raise it - smoothly - to a level that is comfortable for the sonographers to do their job. Because the table is so easy to get on and raise, it helps us when doing stress echoes in that we can get the patient off the treadmill, on the table and into position for imaging in a very short time.

We also like that the Biodex Echo Pro Table has Trendelenburg and Reverse Trendelenburg capability, plus an infinitely adjustable Fowler back. It is comfortable for our patients, looks great in our facility and, thanks to the drop-away sonographer cutout, allows our technicians to get up comfortably close to the patient without bending way over. Because the table is so accommodating, it reduces the time needed to get a quality image.

We highly recommend that other facilities consider the Biodex Echo Pro Table. It has all the features you'll need and it is backed by the Biodex reputation for quality and craftsmanship.”

*- Bonnie Albrecht, PT, Surburban Cardiology, Williamsville, NY*

*“Overall, I'd call the Biodex Ultra Pro Ultrasound {Table} the Cadillac of imaging tables.”*

### Ultra Pro Ultrasound Table



“There are many aspects I like about the Biodex Ultra Pro Ultrasound Table. For starters, the size of this ergonomically designed table is good for patient comfort. It's long enough and wide enough that patients don't end up with their arms or feet hanging off the table.


“The remote control features are another plus. They eliminate the need to physically move, adjust or lift the patient's head or torso; up and down, lateral tilt, Trendelenburg and reverse Trendelenburg - with a 500-pound weight limit, this table can handle and position just about any patient we encounter.

“I also like the arm extension board. During vascular exams, it allows for optimal scanning positioning without the patient's hand draped across your lap. I also find the reverse Trendelenburg positioning for venous exams an asset when a patient is unable to stand.

“Lastly, the locking wheels and floor lock mechanism provide double protection to ensure patient safety.

“Overall, I'd call the Biodex Ultra Pro Ultrasound the Cadillac of imaging tables. It's great for the patient while also lessening the daily wear and tear on my own body.”

*- Lisa LeBlanc, BS, RDMS, RVT, RT(R), (CT),  
Alliance Medical Group, Watertown, CT*





Improving the world's best selling Ultrasound Tables wasn't easy...



...but Biodex was up for the challenge.

Introducing the **All New** line of Ultrasound Tables from Biodex

• ULTRASOUND • ECHOCARDIOGRAPHY • VASCULAR



*Closer patient access*



*Industry's largest cardiac scanning cutout*



*Infinitely adjustable to 80 degrees*



*Maintains table stability*



*Table adjusts down to 23"*

**BIODEX**



# BIODEX

Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9241, Email: [info@biodex.com](mailto:info@biodex.com), [www.biodex.com](http://www.biodex.com)

