



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

QA4 POWERED SURGICAL STRETCHER

Part Number: 465-0200-0000

PRODUCT OVERVIEW

With four remotely-activated powered functions, including a top slide platform, the QA4 Powered Day Surgery Stretchers offer full C-arm coverage and superior surgical access together with a 550 lbs. patient weight limit and exceptional stability.

TECHNICAL SPECIFICATIONS

✓ Length.....	82 ³ / ₄ "
✓ Width.....	25 ³ / ₄ "
✓ Height.....	23 ³ / ₄ " - 39 ¹ / ₂ "
✓ Weight of stretcher.....	350 lbs
✓ Maximum permissible load.....	550 lbs
✓ Trendelenburg/Reverse Trendelenburg.....	20°/12°
✓ Backrest Articulation.....	0-65°
✓ Stretcher Width Stretcher Sides Up.....	30 ¹ / ₄ "
✓ Head Section Articulation.....	+30°/-30°
✓ Lateral Tilt.....	+12°/12°
✓ Head End Top Slide.....	9"
✓ Foot End Top Slide.....	9"

