

# QA4

Surgery Stretcher System

**The complete multi-purpose  
mobile operating system.**



**MEDITEK**

Here for L.I.F.E.

1.800.567.8400

# QA4

Powered Functions

# QA4

Manual Functions

[www.meditek.ca](http://www.meditek.ca)

The QA4 Surgery System, both a patient stretcher and operating table, is a total solution for modern day treatment centres.



- Electrically Powered Functions
- Superior Surgical Access
- Top Slide Platform
- 100% C-Arm Access
- Lateral Tilt
- 550 lbs. Weight Capacity



- Mechanically Powered Functions
- Increased Surgical Access
- Extremely Lightweight & Maneuverable
- 75% C-Arm Access
- Low Cost Solution



Choose from the powered or manual version to suit your needs.

# QA4

Powered Functions

## STANDARD FEATURES


# QA4

Manual Functions

With four remotely-activated powered functions, including a top slide platform, the QA4 Powered stretcher offers full C-arm coverage and superior surgical access together with a 550 lbs. patient weight limit and exceptional stability.



**TRANSPORT & RECOVERY**  
 Stretcher sides for patient comfort and safety  
 Unique retractable steering mechanism engaged for maximum maneuverability



**OPHTHALMICS**  
 Patient platform traversed to head end for access  
 Ophthalmic head rest section



**SHOULDER ARTHROSCOPY**  
 Patient platform traversed to foot end for access  
 Optional arthroscopy backrest

**GYNECOLOGY**  
 Patient platform traversed to foot end for access  
 Leg section removed for access  
 Head section removed for anesthetist access  
 Optional powered leg supports



# QA4

## QA4 Accessory Cart

The complete storage solution for general operating room and QA4 accessories

**GYNECOLOGY**  
 Lightweight leg section removed for access  
 Head section removed for anesthetist access  
 Fitted with optional powered leg supports



**ENT/DENTAL**  
 Patient positioned on standard articulating head section



Excellent maneuverability and easy-to-use hand operated functions make the QA4 Manual stretcher a cost effective option: a patient transport stretcher with the features of an operating table.



**TRANSPORT & RECOVERY**  
 Stretcher sides for patient comfort and safety  
 Unique retractable steering mechanism engaged for maximum maneuverability

1.800.567.8400

www.meditek.ca

**MEDITEK**  
 Here for L.I.F.E.

# QA4 Surgery Stretcher System

We're here for life

STRETCHER SPECIFICATIONS	QA4 POWERED	QA4 MANUAL
	Part No. 21300	Part No. 21310
Powered Functions	✓	
Lateral Tilt	✓	
Head End Top Slide	✓	
Foot End Top Slide	✓	
K8 Pressure Care Mattress	✓	✓
Removable Ultra Lightweight Leg Section	✓	✓
Removable and Rotating Stretcher Sides	✓	✓
Push Button Release and Self Locking Controls	✓	✓
5th Wheel Steering Mechanism	✓	✓
Advanced 4 Caster Braking System	✓	✓
Oxygen Cylinder Trough (Side Loading)	✓	✓
I.V. Pole with Mounting Position	✓	✓
6" Casters	✓	✓
Height (MIN)	23.75"	22.75"
Height (MAX)	39.5"	36.75"
Stretcher Length	82.75"	82.75"
Stretcher Width (With Side Rails)	25.75"	25.75"
Stretcher Width (With Sides)	30.25"	30.25"
Stretcher Width (Brakes Off)	32.75"	32.75"
Mattress Length	80.25"	80.25"
Mattress Width	23.5"	23.5"
Mattress Depth	3"	3"
Base Length	46.75"	41.5"
Trendelenberg	20°	20°
Reverse Trendelenberg	12°	-
Lateral Tilt	+12°/-12°	-
Backrest Articulation	0°-65°	0°-80°
Head Section Articulation	+30°/-30°	+30°/-30°
Head End Top Slide	9"	-
Foot End Top Slide	9"	-
C-Arm Access	100%	75%
Maximum Patient Weight Capacity	550lbs.	350lbs.
Stretcher Weight	350lbs.	275lbs.
Lightweight Leg Section Weight	13lbs.	13lbs.

## Part No. Optional Extras

21312	QA4p Spare I.V. Pole
21313	QA4m Spare I.V. Pole
21315	Sculpted Mattress for Head Support
21320	Ophthalmic Head Support
21340	Removable Shoulder Arthroscopy Backrest
21350	Padded Side Covers - Grey
21351	Padded Sides Covers - Child Print
21355	Articulating Removable Leg Section
21357	Removable Leg Section c/w Full Length Side Rail
21360	K8 Pressure Mattress Complete
21370	QA4 Accessory Cart
21371	QA4 Accessory Cart Spare Trough
21380	Spare Powered Stretcher Hand Control
21385	Spare Powered Stretcher Charging Cable
53556	CD Size Cylinder Support Bracket



Here for L.I.F.E.

1810 Ellice Avenue  
Winnipeg, MB R3H 0B7  
1.800.567.8400  
info@meditek.ca  
www.meditek.ca