

# *QA5 Bariatric Stretcher*




**MEDITEK**

Here for L.I.F.E.

1.800.567.8400 | [info@meditek.ca](mailto:info@meditek.ca) | [www.meditek.ca](http://www.meditek.ca)

# QA5



## Optional Items

- ✓ I.V. Poles
- ✓ Oxygen Tank Holder (Fits into I.V. Socket)

# Bariatric Stretcher

## Technical Specifications

- ✓ Patient weight capacity: 1,000 lbs
- ✓ Overall dimensions: 39" x 84"
- ✓ Patient platform: 37" x 80"
- ✓ Stretcher height range: 20" - 28"
- ✓ Head section angle: 55°
- ✓ Knee angle: 50°
- ✓ Trendelenburg angle: 0 - 10°
- ✓ Reverse trendelenburg angle: 0 - 10°
- ✓ Siderail length: 62"
- ✓ Siderail height above stretcher top: 12"
- ✓ 2 push bars at head and foot ends in addition to fold away push bars at head end
- ✓ Caster size: 8"
- ✓ Floor to base clearance: 9"
- ✓ Floor to lift column clearance: 2"
- ✓ Stretcher Weight: 450 lbs

# Leadership. Innovation. Focus. Excellence.



Cat. No. 11-915