



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

SKYTRON ST19

Part Number: 435-0300-0242

PRODUCT OVERVIEW

SKYTRON's Stellar Series Surgical Lighting provides advanced fixed focus optics within exceptionally durable, lightweight polymer lightheads, creating the perfect balance of performance, maneuverability and dependability. The Vertical Segmented Reflector optics provide shadow controlled illumination and a generous field size. Typical applications include ASC ORs, Labor & Delivery, Exam Treatment and Emergency rooms.

TECHNICAL SPECIFICATIONS

✓ Max Intensity.....	96,800 lux (9,000 fc)
✓ Beam Temperature.....	Less than 25,000 microwatts per sq. centimeter
✓ Color Temperature.....	4000K
✓ Bulb Type.....	Halogen 24V/50W
✓ Fixture Weight.....	72 lbs
✓ Moment Load.....	246 ft-lbs
✓ Electrical.....	120v, 60 Hz, 1 phase

