



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

## **SKYTRON STELLAR ST1919**

Part Number: 435-0300-0243

### **PRODUCT OVERVIEW**

The Skytron Stellar ST1919 utilizes advanced, fixed focus optics within an exceptionally durable lightweight polymer light head, creating the perfect balance of performance maneuverability and dependability. The vertical segmented reflector optics provides shadow free illumination and a generous field size. Its two 19" diameter light heads provide an adjustable intensity of up to 193,600 lux/18,000 foot candles. The light comes with two options; a built in intensity control on the side of the lighthead, or with a wall mounted intensity control (WC).

