



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

## VISTA EXAM LIGHT

Part Number: 435-0300-0254

### TECHNICAL SPECIFICATIONS

- ✓ Large, flood pattern - approximately 6" (150mm) diameter
- ✓ 4,000-ft. candles (43,000 lux) of homogenous light output with over 4,300° Kelvin color temperature at 24" (610 mm) distance
- ✓ Easily replaceable, low cost, long life 12-volt, 75-watt halogen bulb
- ✓ 5" (127mm) glass dichroic reflector
- ✓ IR absorbing, color-correcting glass filter
- ✓ 28" (711 mm) articulating arm with integral compression spring
- ✓ 23" Extension Arm
- ✓ Arm rotates 40° up and down from horizontal plan
- ✓ Lamphead and arm rotate 540 degrees
- ✓ Positioning or optional sterilizable handle mounted on perimeter of lamphead
- ✓ Available in wall, portable floor, single or dual ceiling mounts

