



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

## **MEDITEK QA5 POWERED BARIATRIC STRETCHER**

Part Number: 465-0300-0037

### **PRODUCT OVERVIEW**

Meditek's QA5 Powered Bariatric Stretcher is up to the challenge! Constructed entirely of steel, the QA5 Bariatric Stretcher boasts a 1,000 lbs. weight capacity, as well as heavy duty 8 inch castors. With a built-in CPR release mechanism at the head end, patient left and right side access, the QA5 Bariatric stretcher is geared towards ease-of-use by medical staff and patient comfort.

### **TECHNICAL SPECIFICATIONS**

- ✓ Weight capacity..... 1000 lbs
- ✓ Overall width..... 39"
- ✓ Overall length..... 84"
- ✓ Stretcher top..... 37 x 80"
- ✓ Overall dimensions..... 39" x 83"
- ✓ Stretcher top height..... Adjusts from 20-28"
- ✓ Head section angle..... 55°
- ✓ Knee angle..... 50°
- ✓ Trendelenburg angle..... 0-10°
- ✓ Reverse trendelenburg angle..... 0-10°
- ✓ Siderail length..... 62"
- ✓ Siderail height above stretcher top..... 12"
- ✓ Caster size..... 8"
- ✓ Floor to base clearance..... 9"
- ✓ Floor to lift column clearance..... 2"
- ✓ Stretcher Weight..... 450 lbs

