

## **MEDITEK TABLE SIDE EXTENSIONS**

Part Number: 470-0600-0024; 470-0600-0025; 470-0600-0026; 470-0600-0027

### **PRODUCT OVERVIEW**

Meditek surgical table side extensions are an add-on accessory that allows a general purpose surgical table to be widened to accommodate larger patients. This quickly and easily increases the support surface for your patients and can be attached to the back, seat or leg sections of the table. It may also be used as a support for the arms.

An integral side rail also allows you the ability to attach other accessory attachments you may have on to it, such as armboards or stirrups. Integral clamps ensure that it will be locked securely to the table at all times. Also comes complete with a standard 2" Cushion.

These side extensions can be used on tables that will support patients up to 1000 pounds. Do not exceed the manufacturer listed weight capacity for the surgical table.

- ✓ 470-0600-0024 Table Side Extension – 4" x 20"
- ✓ 470-0600-0025 Table Side Extension – 4" x 15"
- ✓ 470-0600-0026 Table Side Extension – 4" x 12"
- ✓ 470-0600-0027 Table Side Extension – 4" x 23"

