



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

## **MEDITEK AS-4000 ARMBOARD**

Part Number: 470-1310-0019

### **PRODUCT OVERVIEW**

- ✓ Premium quality standard armboard with a phenolic operating surface
- ✓ Armboard is 26" long
- ✓ Rotates 180° and can be positioned in 21 locking increments
- ✓ Easily attaches to all North American surgical tables
- ✓ A lightweight and strong alloy is used for the gear, eliminating the possibility for broken gear teeth
- ✓ Cushions are available in 1", 2", 3" or 4" sizes

