



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

MEDITEK AS-6000 ARMBOARD

Part Number: 470-1310-0029

PRODUCT OVERVIEW

- ✓ 6"W x 20"L Phenolic operating surface
- ✓ Permits complete range of positioning flexibility
- ✓ Stainless steel support bar
- ✓ Comes with 2" Pressure Management pad
- ✓ Provides arm support for Lateral, Prone, Neuro,
- ✓ Requires clark socket or super clamp for mounting to siderail of surgical table - Sold separately

