



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

MEDITEK LEVITATOR STIRRUPS

Part Number: 470-1330-0059

PRODUCT OVERVIEW

The Meditek Levitator Stirrups are the ultimate in patient positioning, for all types of gynecology and urology requirements. The blade type design means that they will never slip out of place after they have automatically locked into position. A gas shock assists you in effortlessly positioning your patient quickly and efficiently.

The Meditek Levitator Stirrups are single handedly adjusted and positioned from outside the sterile field. Patient comfort is guaranteed with our memory foam gel pad that fully encapsulates the foot. The standard weight capacity of this pair is 600 pounds.

