

## **BIODEX SURGICAL C-ARM TABLE - 840**

Part Number: 475-0100-0003

### **PRODUCT OVERVIEW**

Featuring a large radiolucent area, motorized actuation of height, lateral roll and Trendelenburg motions, the Pain Management C-Arm Table is ideal for pain care applications.

### **TECHNICAL SPECIFICATIONS**

✓ Length.....	104" (260 cm)
✓ Width.....	26" (66 cm)
✓ Maximum permissible load (Lift).....	500 lb (227 kg)
✓ Radiolucent Area.....	51" l x 24" w (128 x 61 cm)
✓ Adjustments:	
✓ Height.....	33.5" - 43.5" (84 - 109 cm)
✓ Tabletop X Motion (head-to-toe float).....	20" (51 cm)
✓ Tabletop Y Motion (side-to-side float).....	10" (25.4 cm)
✓ Trendelenburg.....	20°
✓ Reverse Trendelenburg.....	20°
✓ Lateral Roll.....	0°/20°

