

Steris SurgiGraphic 6000 Imaging Surgical Table

Part Number: 475-0100-0009

PRODUCT OVERVIEW

The SurgiGraphic surgical table enables radiology quality imaging for minimally-invasive vascular, orthopedic, or other procedures requiring fluoroscopy; lateral and longitudinal float, Trendelenburg and reverse Trendelenburg, lateral tilt, height adjustment, and the most useful combinations of these movements.

TECHNICAL SPECIFICATIONS

- ✓ Length.....92" (2340 mm)
- ✓ Width.....21" (530 mm)
- ✓ Maximum permissible load (Lift).....500 lbs 227 kg
- ✓ Maximum permissible load (Articulation).....500 lbs 227 kg
- ✓ Weight of Table.....900 lbs 449 kg
- ✓ Adjustments:
 - ✓ Height.....32" - 42" (812 - 1067 mm)
 - ✓ Lateral tilt (left/right).....16°
 - ✓ Trendelenburg.....16°
 - ✓ Reverse Trendelenburg.....16°

