



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

MEDSTONE ELITE TE

Part Number: 475-0100-0018

PRODUCT OVERVIEW

The Medstone TM-5 Imaging Surgical Table is designed for high patient throughput, utilizing the highest quality components and designed for ease of use by the medical staff, with patient comfort in mind.

The table top has been designed to give maximum unobstructed Imaging, allowing 360° rotation radiation free coverage around the longitudinal axis of the carbon fibre table-top which makes it ideal for 3D imaging studies.

TECHNICAL SPECIFICATIONS

- ✓ Weight Capacity.....500 lbs
- ✓ Surface Dimensions.....22" x 82"
- ✓ Unobstructed Tabletop Imaging Area.....58" (1473mm)
- ✓ Adjustments:
 - ✓ Height.....31" - 43"
 - ✓ Lateral tilt (left/right).....15°
 - ✓ Trendelenburg.....15°
 - ✓ Reverse Trendelenburg.....15°

