



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

MAQUET 1130

Part Number: 475-0200-0000; 475-0200-0001

PRODUCT OVERVIEW

The Maquet 1130 surgical table is ideal for radiography, intraoperative image intensification and is C-Arm compatible. It is mobile with a battery powered electrically controlled hydraulic drive, and the tabletop is divided into up to 8 sections. The padded plates for the upper/lower sections and seat plate can be individually removed. The operating table has a T-shaped base with large, easy-rolling, ball-bearing castors for lateral and transverse motion.

TECHNICAL SPECIFICATIONS

- ✓ Weight Capacity.....300 lbs
- ✓ Adjustments:
 - ✓ Height.....28" - 43"
 - ✓ Lateral tilt (left/right).....25°
 - ✓ Trendelenburg.....30°
 - ✓ Reverse Trendelenburg.....30°
 - ✓ Leg up/down.....25°/-75°
 - ✓ Back up/down.....65°/45°
 - ✓ Flex.....Adjustable
 - ✓ Reflex.....+10 back/+15 seat
 - ✓ Return to level.....Yes

