

MAQUET ALPHASTAR 1132

Part Number: 475-0200-0002

PRODUCT OVERVIEW

The Alphastar 1132 Surgical Table has all of the essential functions to be truly versatile for many surgical specialties.

TECHNICAL SPECIFICATIONS

- ✓ Length..... 71.3" (1810 mm)
- ✓ Length (with headrest)..... 82.7" (2100 mm)
- ✓ Width (back/seat plate)..... 20.5/21.2" (520/535 mm)
- ✓ Width (with side rails) (back/seat plate)..... 22/23" (560/585 mm)
- ✓ Weight of operating table..... 715 lbs (325 kg)
- ✓ Maximum permissible load..... 500 lbs (225 kg)
 - ✓ (1132.03) Normal..... 800 lbs (360 kg)
 - ✓ (1132.03) With upper back extension plate..... 500 lbs (225 kg)
- ✓ Adjustments via hand control:
 - ✓ Height (including cushions)..... 28.9 - 46.1" (735 - 1170 mm)
 - ✓ Lateral tilt (left/right)..... 18°
 - ✓ Trendelenburg..... 30°
 - ✓ Reverse Trendelenburg..... 30°
 - ✓ Back up/down..... 68°/25°
 - ✓ Leg up/down..... 80°/105°
 - ✓ Flex/Reflex
 - ✓ Base braking (lock/unlock)
 - ✓ Patient orientation (normal/reverse)
 - ✓ Return to level (anatomically correct return to horizontal position from articulated positions)

