

MAQUET ALPHAMAXX 1133

Part Number: 475-0200-0003

PRODUCT OVERVIEW

The Maquet 1133.02 Alphamaxx is a versatile top slide surgical table designed to enhance patient safety and comfort. It can withstand extreme loads in every configuration and can be set up for any surgical discipline in just a few easy steps.

TECHNICAL SPECIFICATIONS

- ✓ Length..... 86" (2180 mm)
- ✓ Width..... 21" (533 mm)
- ✓ Weight of operating table..... 693 lbs (315 kg)
- ✓ Maximum permissible load (Lift & Articulation)..... 800 lbs (363kg)
- ✓ Maximum permissible load (Side Tilt)..... 800 lbs (363kg)
- ✓ Adjustments via hand control:
 - ✓ Height..... 24 - 41" (610 - 1041 mm)
 - ✓ Lateral tilt (left/right)..... 20°
 - ✓ Trendelenburg..... 30°
 - ✓ Reverse Trendelenburg..... 30°
 - ✓ Back up/down..... 80°/40°
 - ✓ Leg up/down..... 10°/90°
 - ✓ Flex/Reflex..... 20°/55°
 - ✓ Slide..... 9" (229mm)
- ✓ Base braking (lock/unlock)
- ✓ Return to level (anatomically correct return to horizontal position from articulated positions)

