

SKYTRON 1602 ESSENTIA SURGICAL TABLE

Part Number: 475-0200-0006

PRODUCT OVERVIEW

The 1602 Essentia allows the table top to rotate independent of the base, allowing comfortable surgical access. Featuring versatile positioning options such as 90° back-up, this table is a great value for operating rooms and surgery centers.

TECHNICAL SPECIFICATIONS

- ✓ Weight of operating table.....500 lb (227kg)
- ✓ Maximum permissible load.....500 lb (227kg)
- ✓ Adjustments via hand control:
 - ✓ Height.....27 ³/₁₆ to 40 ¹⁵/₁₆" (690 - 1,040 mm)
 - ✓ Lateral tilt (left/right).....23°
 - ✓ Trendelenburg/Reverse.....28°
 - ✓ Back up/down.....90°/30°
 - ✓ Leg down.....100°
 - ✓ Return to level (anatomically correct return to horizontal position from articulated positions)
- ✓ Manual adjustments:
 - ✓ Top Rotation.....180°
 - ✓ Head Section.....60° up/90° down

