



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

## SKYTRON 3602 ULTRASLIDE

Part Number: 475-0200-0010

### PRODUCT OVERVIEW

Skytron's UltraSlide 3602 table is uniquely designed to address each of your OR's needs. Its industry-leading 23 inch top slide allows you to easily move patients through the imaging window and pre-programmed positioning options allow for quick table set-up at any point during the procedure. Full body imaging can be accomplished without making risky patient maneuvers, making the UltraSlide 3602 a safe and easy to use solution. The 1000 lb. lift and 800 lb. articulation capacity is never compromised by top slide, meaning that your OR can more efficiently support a greater variety of procedures.

### TECHNICAL SPECIFICATIONS

- ✓ Length w/o head rest.....69" (1750 mm)
- ✓ Length with head rest.....82" (2080 mm)
- ✓ Width.....19½" (500 mm)
- ✓ Width with side rails.....22" (560 mm)
- ✓ Weight of operating table.....800 lbs (365 kg)
- ✓ Maximum permissible load (Lift).....1000 lbs (455 kg)
- ✓ Maximum permissible load (Articulation).....800 lbs (365 kg)
- ✓ Adjustments:
  - ✓ Height.....24-42½ (610-1080 mm)
  - ✓ Lateral tilt (left/right).....30°
  - ✓ Trendelenburg.....30°
  - ✓ Reverse Trendelenburg.....30°
  - ✓ Power Kidney Bridge.....5¾ (145 mm)
  - ✓ Back up/down.....90°/40°
  - ✓ Leg down.....15°/100°
  - ✓ Return to level.....Yes
  - ✓ Top slide/head.....9½" (240 mm)
  - ✓ Top slide/foot.....13½" (340 mm)
  - ✓ Flex/Reflex.....40°/90°
  - ✓ Base braking (lock/unlock).....Yes
  - ✓ Beach Chair.....Yes

