



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

AMSCO 2080 RC

Part Number: 475-0200-0022

TECHNICAL SPECIFICATIONS

✓ Length (with headrest).....	84" (2117 mm)
✓ Width.....	20" (508 mm)
✓ Weight of operating table.....	690 lbs (313 kg)
✓ Maximum permissible load (Lift).....	300 lbs (136 kg)
✓ Maximum permissible load (Articulation).....	300 lbs (136 kg)
✓ Adjustments:	
✓ Height.....	27 - 45" (686 - 1134 mm)
✓ Lateral tilt (left/right).....	10°
✓ Trendelenburg.....	20°
✓ Reverse Trendelenburg.....	20°
✓ Kidney Bridge.....	4½" (113 mm)
✓ Back up/down.....	90°/10°
✓ Leg up/down.....	0°/105°
✓ Flex.....	15° Back/20° Seat
✓ Reflex.....	10° Back/15° Seat

