



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

STERIS 3085 SP

Part Number: 475-0200-0025

PRODUCT OVERVIEW

The 3085 SP Surgical Table provides complete flexibility for patient positioning, no matter how difficult the procedural or imaging requirements. Its proven performance ensures increased productivity, ease of use and long term reliability. The 3085 SP handles virtually all the surgical specialties you expect a general surgical table to handle, including urology and neurology, pediatrics and bariatrics, ENT and endovascular, and even orthopedics, including skeletal traction, with accessories.

TECHNICAL SPECIFICATIONS

- ✓ Length.....76" (1930 mm)
- ✓ Width.....20" (508 mm)
- ✓ Weight of operating table.....737 lbs (334 kg)
- ✓ Maximum permissible load (Lift & Articulation).....1000 lbs (454 kg)
- ✓ Maximum permissible load (Side Tilt).....500 lbs (227 kg)
- ✓ Adjustments via hand control:
 - ✓ Height.....27 - 44" (686 - 1118 mm)
 - ✓ Lateral tilt (left/right).....18°
 - ✓ Trendelenburg.....25°
 - ✓ Reverse Trendelenburg.....25°
 - ✓ Manual Kidney Bridge.....4" (150 mm)
 - ✓ Back up/down.....55°/25°
 - ✓ Leg up/down.....80°/105°
 - ✓ Flex/Reflex.....20°/25°
 - ✓ Base braking (lock/unlock)
 - ✓ Return to level (anatomically correct return to horizontal position from articulated positions)

