



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

## **SKYTRON 3600B ULTRASLIDE**

Part Number: 475-0200-0028

### **PRODUCT OVERVIEW**

The 3600B ULTRASLIDE Surgical Table provides 23 in. Top Slide, full table features and optimal patient positioning for C-arm access, including powerful 1,000 lb. lift and 800 lb. articulation. Removable back and leg sections and Auto Beach Chair position provide flexible patient positioning options for every procedure.

### **TECHNICAL SPECIFICATIONS**

- ✓ Length..... 69" (1750 mm)
- ✓ Length with head rest..... 82" (2080 mm)
- ✓ Width..... 19.5" (500 mm)
- ✓ Width (with side rails)..... 22" (560 mm)
- ✓ Weight of operating table..... 800 lbs (365 kg)
- ✓ Maximum permissible load (Lift)..... 1000 lbs (455 kg)
- ✓ Maximum permissible load (Articulation)..... 800 lbs (365 kg)
- ✓ Adjustments via hand control:
  - ✓ Height..... 24 - 42.5" (610 - 1080 mm)
  - ✓ Lateral tilt (left/right)..... 30°
  - ✓ Trendelenburg..... 30°
  - ✓ Reverse Trendelenburg..... 30°
  - ✓ Power Kidney Bridge..... 5.75" (145 mm)
  - ✓ Back up/down..... 90°/40°
  - ✓ Leg up/down..... 15°/100°
  - ✓ Top slide/head..... 9.5" (240 mm)
  - ✓ Top slide/foot..... 13.5" (340 mm)
  - ✓ Flex/Reflex..... 40°/90°
- ✓ Base braking (lock/unlock)
- ✓ Beach Chair
- ✓ Return to level (anatomically correct return to horizontal position from articulated positions)

