

SHAMPAINE 4900B

Part Number: 475-0200-0029

PRODUCT OVERVIEW

The Shampaine 4900B is a remotely controlled, full functional surgical table that provides flexible articulation for patient positioning. A longitudinal slide feature provides full-length image amplification and a cassette tunnel is incorporated into each of the five table sections. The 4900B is particularly adaptable for procedures requiring mobile image amplification, eliminating the need for additional awkward accessories. No other table provides greater C-arm accessibility for laparoscopic/endoscopic surgery. The 4900B is battery operated.

TECHNICAL SPECIFICATIONS

✓ Length.....	71.5" (1842 mm); 88" with removable leg section
✓ Width.....	19.75" (502 mm); 21.75" with side rails
✓ Height.....	28 - 43" (711 - 1092 mm)
✓ Lateral tilt (left/right).....	18°
✓ Trendelenburg.....	25°
✓ Reverse Trendelenburg.....	25°
✓ Leg up/down.....	70°/95°
✓ Back Section.....	30°/90°
✓ Kidney Elevator.....	Extendable up to 2.5" (64 mm)
✓ Headrest Section.....	45°/90° (in 15° increments)
✓ Battery Operation	

