



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

SKYTRON 6500B

Part Number: 475-0200-0035

PRODUCT OVERVIEW

The 6500B Surgical Table was designed for flexible patient positioning, unobstructed Image Amplification and enhanced surgeon comfort. Its uncommon flexibility makes the 6500B unequalled for general purpose surgery. The exclusive 180° rotating top feature allows unobstructed C-arm and microscope Coverage. Additionally, it eliminates conflicts between surgeon and anesthesiologist during head-end procedures, and provides ample leg room for surgeons who prefer to sit during procedures.

TECHNICAL SPECIFICATIONS

- ✓ Length..... 65.5" (1665 mm)
- ✓ Length (with head rest)..... 79" (2010 mm)
- ✓ Width..... 20" (504 mm)
- ✓ Width (with side rails)..... 21.75" (548 mm)
- ✓ Max. permissible load (Lift)..... 500 lbs. (227 kg)
- ✓ Max. permissible load (Articulation)..... 500 lbs. (227 kg)
- ✓ Top rotation..... 180°
- ✓ Adjustments via hand control:
 - ✓ Height..... 26 - 47" (655 - 1184 mm)
 - ✓ Lateral tilt (left/right)..... 30°
 - ✓ Trendelenburg..... 30°
 - ✓ Reverse Trendelenburg..... 30°
 - ✓ Power Kidney Bridge..... 5" (126 mm)
 - ✓ Back up/down..... 90°/40°
 - ✓ Leg up/down..... 20°/105°
 - ✓ Flex/Reflex..... 40°/90°
 - ✓ Base braking (lock/unlock)
 - ✓ Return to level (anatomically correct return to horizontal position from articulated positions)
- ✓ Battery powered

