

AMSCO 2080 L

Part Number: 475-0200-0048

PRODUCT OVERVIEW

Amsco 2080 General Surgical Tables are built with quality and are built to last. The Amsco 2080 General Surgical Tables are easy to use, versatile and well within your budget.

The Amsco 2080 General Surgical Tables are designed to accommodate a maximum 300 pound load, distributed as follows:

- ✓ Head section 25 lbs.
- ✓ Back section 100 lbs.
- ✓ Seat section 140 lbs.
- ✓ Foot section 35 lbs.

TECHNICAL SPECIFICATIONS

- ✓ Adjustments:
 - ✓ Height..... 27 - 45"
 - ✓ Lateral tilt (left/right)..... 10°
 - ✓ Trendelenburg..... 20°
 - ✓ Reverse Trendelenburg..... 20°
 - ✓ Kidney Up..... 4½"
 - ✓ Back up/down..... 90°/-10°
 - ✓ Flex..... -15 back/-20 seat
 - ✓ Reflex..... +10 back/+15 seat
 - ✓ Return to level..... No

