



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

## STERIS 5085

Part Number: 475-0200-0054

### PRODUCT OVERVIEW

The 5085 Surgical Table can adapt to virtually every surgical specialty. It boasts an exceptionally low height capability, high weight limit and large imaging area. The pre-programmed hand control allows for quick, one touch patient positioning for beach chair or rotation.

The 5085 can easily be upgraded to the 5085 SRT, which features drive ability so the table can also transport patients. It also allows 360° rotation during surgery with no support posts or weight calculations.

### TECHNICAL SPECIFICATIONS

- ✓ Length..... 81" (2057 mm)
- ✓ Width..... 20" (508 mm) (28° with extenders)
- ✓ Weight Capacities
  - ✓ Patient support (Raise & Lower)..... 1,200 lbs. (544 kg)
  - ✓ Maximum permissible load (Lift & Articulation)..... 1,000 lbs (454 kg)
  - ✓ Maximum permissible load (slide)..... 600 lbs (272 kg)
- ✓ Table Rotation..... 360°
- ✓ Weight of operating table..... 820 lbs. (372 kg)
- ✓ Head Section..... +90° to -90° (15° increments)
- ✓ Adjustments via hand control:
  - ✓ Height..... 22.8 – 42.5" (581 - 1080 mm)
  - ✓ Lateral tilt (left/right)..... 20°
  - ✓ Trendelenburg/Reverse..... 30°
  - ✓ Top Slide..... 9" (227 mm) to head, 9" to foot
  - ✓ Powered Kidney Bridge..... 4" (150 mm)
  - ✓ Back up/down..... 80°/40°
  - ✓ Leg up/down..... 0°/105°
  - ✓ Flex/Reflex..... 140°/100°
  - ✓ Base braking (lock/unlock)
  - ✓ Return to level (anatomically correct return to horizontal position from articulated positions)

