

Diamond Series

Ceiling Column



Medical Design's retractable or fixed Ceiling Column offers an ergonomic solution for the distribution of medical gases, electrical services and communication devices in operating room suites and emergency wards. This concept allows for a better management of workspace in the patient care area.

The DIAMOND series makes it possible to optimize the quantity of services while leaving a maximum of accessibility to the user. Our units are offered in stainless steel or in epoxy coated steel with a light almond color finish.

Diamond Series

Ceiling Column

Epoxy coated finish or stainless steel

Epoxy coated finish or stainless shroud

Rear access panel

Electricity, data and communication services are available on the side for a maximum quantity of 18

Medical Gas services are available under the unit up to a quantity of 10 and located as per customers needs

Applications

Emergency Wards

Medical Imagery

Operating Rooms

Traumatology

Anesthesia

Endoscopy

Product Benefits

- These columns are offered in the electric retractable version or fixed model.
- The Diamond is ergonomically designed to accommodate up to 28 services.
- Available with integrated CO2 or nitrogen control panel.
- The unit allows installation of a peripheral rail under the column to add related accessories.
- When raised, the electric retractable column allows more space. It can easily be lowered to give trouble-free hook-up to services.
- Electrical devices are pre-wired and medical gas outlets are terminated to a single point connection at the top of the unit.

