

Epidural Positioning Devices



M E D I T E K

Here for L.I.F.E.

1.800.567.8400 | info@meditek.ca | www.meditek.ca

Manual & Electric EPDs

The Epidural Positioning Device commonly known as the EPD, is a revolutionary device for positioning patients safely and securely during epidural procedures.

- Reduces strain and injury risk to staff
- Correctly positions the patient by cervical, thoracic & lumbar flexion, maintaining a solid and stable position
- Provides patient comfort with less invasion of their personal space compared to traditional positioning methods

	Manual	Electric
EPD Specifications		
Weight	12 lbs	110 lbs
Load Rating	400 lbs	400 lbs
Arm Rest Range	32" - 43"	6" Range
Arm Rest Size	16.5" x 10"	16.75" X 4.5"
Chest Piece Range	5.5" - 10"	6" Range
Chest Piece Depth	3"	N/A
Chest Piece Size	8" x 8"	6.5" x 6.5"
Shipping Weight	42 lbs	170 lbs.
Shipping Box	32" x 18" x 15"	36.5"x38.5"x50"
LDS Specifications (base)		
Weight	42 lbs	N/A
Foot Plate Size	20" x 12"	
Foot Plate Range	10.5" x 19.5"	
Width	20.5"	
Length	23"	
Shipping Weight	50 lbs	
Shipping Box	29" x 26" x 21"	



Increase efficiency in operating rooms, labour and delivery suites, pain centres and any other facility/department where epidurals and spinals are performed.

Manual EPD

180° Adjustable Face Rest

Tilts to attain and maintain proper cervical flexion

6 Position Adjustable Arm Rest

Adjustable range from the floor 32" - 43";
Allows patient to rest arms comfortably and provides added stability

Adjustable Torso Support

Increase lumbar and thoracic flexion with simple adjustments; chest support immobilizes back

*Optional flat torso cushion



Positive Locking Clamp

Easily adjusts; locks with one turn

7 Position Adjustable Foot Rest

Adjustable range from the floor 10.5" - 19.5"
accommodates a wide range of body sizes

High Quality Casters

Easily lock into place



Electric EPD

180° Adjustable Face Rest

With repositionable pillow for use during cervical as well as lumbar and thoracic placement

Adjustable Torso Support

Increase lumbar and thoracic flexion with simple adjustments; chest support immobilizes back

Hand Control

Ergonomic hand control with convenient hanger (Controls both midsection and foot rest)

Adjustable Arm Rest

Allows patient to rest arms comfortably, provides added stability

Battery

Rechargeable battery lasts up to 100 cycles (40 lifts and lowerings for the JUMBOsystem with fully charged battery)

Adjustable Foot Rest

Adjusts to the proper height with hand control to accommodate a wide variety of patients

High Quality Casters

Low profile castor locks, conveniently tucked out of the way of pant legs and feet, yet easy to access for locking down and releasing





Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.



- ✓ Cushions & Mattresses
- ✓ Equipment Booms
- ✓ Gel Pads & Positioners
- ✓ Linen Solutions

- ✓ Nurse Documentation Centers
- ✓ OR Products
- ✓ Patient Transfer
- ✓ Room Equipment

- ✓ Sterilization & CSR Products
- ✓ Stretchers & Accessories
- ✓ Surgical & Exam Lights
- ✓ Surgical Tables & Accessories

Cat. No. 02-917