
Providing Quality

Capital Equipment

For Healthcare Professionals



Stretcher and Patient Transport Overview

MEDITEK

Here for L.I.F.E.

Stretcher and Patient Transport Overview

1100 Fixed Height Stretcher

- 500 lb weight capacity
- Fixed height available in 3 options: 28", 30", 32"
- Pneumatic assist head section
- Four corner brake/steer pedals
- Full length tuck away side rails



1100 Fixed Height Stretcher



4500 GYN
Stretcher

4500 GYN Stretcher

- Sliding top deck, converts from stretcher to chair
- Pneumatic assist head and leg sections
- Integrated telescoping heel stirrups that fold away when not in use
- Folding collapsible side rails
- Telescoping IV pole at head end
- Height adjustment of 24" to 37"
- Raise/lower pedals located at head end and on both sides of stretcher

4500 Stretcher

- 500 lb weight capacity
- Hydraulic height adjustment of 22" to 36"
- Optional four-sided raise/lower pedals
- Pneumatic assist head section
- Easy clean base with integrated storage and oxygen tank holder
- Folding collapsible side rails
- Available in narrow 29" wide version



4500 Stretcher

Easy Lift Transfer Stretcher

- 300 lb weight capacity
- Crank raise/lower control
- Height adjustment of 30" to 42"
- Variety of sheet options available for: bathing, cysto/OBS, etc.
- Easy to use self locking buckles



Easy Lift
Transfer Stretcher

M E D I T E K



5000 Stretcher

5000 Stretcher

- 500 lb weight capacity
- Bariatric 650 lb version available
- Hydraulic height adjustment of 22" to 36"
- Optional 4 sided raise/lower pedals
- Pneumatic assist head section
- Easy clean base with integrated storage and oxygen tank holder
- Full length tuck away side rails

5000 Cadaver Stretcher

- 500 lb weight capacity
- Bariatric 650 lb version available
- Hydraulic height adjustment of 22" to 32"
- Head and foot end raise/lower pedals
- Lightweight removable frame with vinyl covered canopy



5000 Cadaver
Stretcher

Surgical Tables
Table Accessories
Surgical Lights
Cushions
Armboards
Stirrups
Gel Pads
Stretchers
Patient Positioners
Restraint Straps
Hand Surgery Tables
Preventative Maintenance
And Much More!



Here for L.I.F.E.

1810 Ellice Avenue
Winnipeg, MB R3H 0B7
1.800.567.8400
info@meditek.ca
www.meditek.ca