

SMART Rollboards

SMART transfer boards remove the need to lift when performing patient transfers. SMART transfer boards are manufactured and designed so that staff do not have to use extra effort, which makes their job easier and safer.

Features and Benefits:

- High performance
- Radiolucent
- Non-toxic materials
- Flame retardant surface
- Anti-bacterially surface
- Safe and comfortable for patients
- No Lifting
- Prevents back injuries to care workers
- Latex-free/skin friendly
- Easy to clean and disinfect
- Disposable covers available
- Tear-resistant nylon cloth
- Bridges gaps up to 20 cm
- Lightweight, easily suspendable
- MRI safe
- 3-year warranty

Easy Patient Transfers

Avoid injuries to patients or staff during transfers.



MEDITEK

Transfer patients in 3 easy steps:

Turn the patient onto the side, slide the board under the body and return the patient to the original position.



Gently push the patient at the shoulder and hip, without applying force, and transfer.



Turn the patient onto the side. Remove the Rollboard starting at the feet.



SMART rollboards are made of cutting-edge materials and manufacturing processes to provide a seamless transition from one surface to another. With our devices, there is no need to lift a patient. Minimal effort is required to perform the transfer. This reduces the strain on the caregiver's body while providing a comfortable experience for the patient.



Product Name	Specifications	Product No.
	<p>Standard SMART Rollboard Suggested Use: Radiology, OR, ICU, General Ward</p>	<p>Size: 177 x 50 cm / 70 x 20 inch foldable</p> <p>Weight: 2.7 kg</p>
	<p>Surgiboard Long Narrow SMART Rollboard Suggested Use: OR and X-Ray</p>	<p>Size: 120 x 40 cm / 48x 16 inch non-fold- able</p> <p>Weight: 1.5 kg</p>
	<p>MRI / CT SMART Rollboard Suggested Use: To or From Stretchers, in CT and MRI</p>	<p>Size: 177 x 40 cm / 70x 16 inch foldable</p> <p>Weight: 2.5 kg</p>
	<p>Mini Short SMART Rollboard Suggested Use: Urol- ogy, OR, Gynecolo- gy, Recovery Rooms, Emergencies, Pediat- ric, and Birthing</p>	<p>Size: 91 x 50 cm / 36x 20 inch non-fold- able</p> <p>Weight: 1.4 kg</p>



Product Name	Specifications	Product No.	
	<p>Surgiboard Short Narrow SMART Rollboard Suggested Use: OR, X-Ray, Ambulance, and Pediatric</p>	<p>Size: 100 x 40 cm / 40x 16 inch non-fold- able</p> <p>Weight: 1.3 kg</p>	470-1400-0049
	<p>Mini Long SMART Rollboard Suggested Use: OR, X-Ray, and Birthing</p>	<p>Size: 100 x 50 cm / 40x 20 inch non-fold- able</p> <p>Weight: 1.5 kg</p>	470-1400-0050
	<p>Surgi Wide SMART Rollboard Suggested Use: OR, X-Ray, and Pediatric</p>	<p>Size: 120 x 50 cm / 48x 20 inch non-fold- able</p> <p>Weight: 1.5 kg</p>	470-1400-0051

