

Surgical Table Accessories



M E D I T E K

Here for L.I.F.E.

Welcome!



Here for L.I.F.E.

Since 1981, Canadian healthcare facilities have been receiving outstanding customer care in their dealings with Meditek. Patients and front line employees alike appreciate the reliability Meditek offers through superior quality healthcare products and services.

You care about your patients' safety and well-being. Our job is to help you facilitate that care to the highest standard, while exceeding your expectations.

Personalized care and support are the cornerstones of what Meditek is. Coast to coast territory sales managers are ready to assist you in any way they can. With every contact, you will receive their full attention, vast medical products expertise and Meditek's commitment to provide you with a solution that fits your needs.

You know your needs better than anyone. By listening to frontline user input, Meditek has become an established designer, manufacturer and supplier of a wide variety of medical products, from mattresses to stretchers.

Our commitment to providing your team with the equipment they need in today's fast changing environment gives **Here for L.I.F.E.** real meaning!

Leadership. Innovation. Focus. Excellence.

The Meditek Team



pg. 7



pg. 12



pg. 4



pg. 8



pg. 17

Table of Contents

Anesthesia Screens	4
Armboards	5
Arm & Hand Surgery Tables	7
Clamps	8
Foot & Side Extensions	9
Leg Supports	10
Positioning Devices	14
Restraint Straps	15
Shoulder Chairs	17
Storage Carts & Holders	18
Urology Systems	19

4 *Anesthesia Screens*



Universal Anesthesia Screen without Clark Socket | 470-0100-0006

Universal Anesthesia Screen | 470-0100-0007

- ✓ Comes with one 470-0300-0028 Clark Socket clamp

Form-A-Screen | 470-0100-0008

- ✓ Flexible design
- ✓ Requires two 470-0300-0028 Clark Socket clamps

Form-A-Screen | 470-0100-0009

- ✓ Flexible design
- ✓ Comes with two 470-0300-0028 Clark Socket clamps



Standard Armboard



- ✓ Rotates 180° and can be positioned in 21 locking increments
- ✓ 26" long phenolic operating surface

Meditek Armboard without Cushion | 470-1310-0019

Meditek Armboard Standard 2" Cushion | 470-0400-0286



Meditek Armboard Standard 3" Cushion | 470-0400-0288

Meditek Armboard Standard 4" Cushion | 470-0400-0289

Arm Support with Secure Lock Bracket without Cushion | 470-1310-0028

- ✓ Height adjustable to accommodate all thicknesses of cushions
- ✓ 26" long phenolic operating surface
- ✓ Requires 470-0400-0286 Cushion
- ✓ No interference during imaging and X-ray procedures
- ✓ Securely locks to operating room table with integrated clamp



Multi-Position Armboard | 470-1310-0029

- ✓ 6"W x 20"L phenolic operating surface
- ✓ Comes with 2" pressure management cushion
- ✓ Provides arm support for either lateral or prone procedures
- ✓ Requires 470-0300-0028 Clark Socket clamp for mounting to side rail of table

6 Armboards



Lateral Armboard | 470-1310-0031

- ✓ Requires 470-0400-0286 Cushion

Arm Support Guard | 470-0920-0137

- ✓ 7.5" x 17" size
- ✓ Works on any surgical table section
- ✓ Ample room for securely holding any outer extremity rigidly in place
- ✓ Can be used on any size of patient
- ✓ Very easy to clean



Arm Support Guard with Relaxa Gel Pad | 470-0920-0088

- ✓ Includes re-usable viscoelastic polymer gel pad ensures optimal patient comfort
- ✓ Reduces the risk of patient injury, including decubitus ulcers, caused by sheer friction



Kraus Arm Support without Clark Socket | 470-1310-0035

- ✓ Stainless steel arm cradle measures 6" x 20" overall
- ✓ Perfect for arm chair, prone, lateral and neuro procedure positioning
- ✓ Can be used as a well-leg positioner

Arm & Hand Surgery Tables 7



Legless Lightweight Arm & Hand Table | 470-0200-0007

- ✓ Complete with 1" cushion
- ✓ Weighs only 7 lbs
- ✓ Size 16" x 30"
- ✓ Requires two 470-0300-0021 side rail clamps

Free Access Arm and Hand Surgery Table | 470-0200-0008

- ✓ Complete with 2" cushion and clamps
- ✓ Size 12" x 30"



Legless Arm & Hand Surgery Table | 470-0200-0015

- ✓ Complete with 2" cushion and clamps
- ✓ Size 18" x 32"

Single Leg Support Arm & Hand Surgery Table | 470-0200-0009

- ✓ Complete with 2" cushion and clamps
- ✓ Size 18" x 32"
- ✓ Single leg support rod for added stability



Hourglass Arm & Hand Surgery Table | 470-0200-0013

- ✓ Complete with 2" cushion and clamps
- ✓ Size 15" x 32"
- ✓ Single leg support rod for added stability

8 *Clamps*



Side Rail/Foot Extension Clamp | 470-0300-0021

- ✓ Accommodates all rectangle bar accessories up to 1/4" x 1"
- ✓ For safe, secure positioning of straight bar accessories
- ✓ Stainless steel for enhanced durability

Premium Side Rail/Foot Extension Clamp | 470-0300-0023

- ✓ Lightweight solid construction aluminum
- ✓ Accepts standard 1" x 3/8" blade style accessories



Meditek Easy Lock Socket Clamp | 470-0300-0025

- ✓ Allows smooth and easy positioning of surgical table accessories
- ✓ 3 Easy pull-to-release lever
- ✓ Lithotomy wheel for accessory angle adjustments
- ✓ Fits up to 1-1/16" diameter round posts

Clark Socket Clamp | 470-0300-0028

- ✓ Allows smooth and easy positioning of surgical table accessories
- ✓ Accommodates all accessories up to 3/4" diameter





Stainless Steel Foot Extensions | 470-0400-0302

- ✓ Can be used for extending your O.R. Table or as a reverse trendelenburg foot support platform
- ✓ Includes 2" cushion and two 470-0300-0021 side rail clamps
- ✓ 8" length

Meditek O.R. Table Side Extension

- ✓ Includes 2" cushion
- ✓ Integral locking clamps
- ✓ Adds 4" of width to any O.R. table
- ✓ Side rail bar

20" length | 470-0600-0024

15" length | 470-0600-0025

12" length | 470-0600-0026

23" length | 470-0600-0027



10 *Leg Supports*



Altima Leg Holders | 470-1330-0056

- ✓ Stainless steel adjustable extension bar for flexible positioning
- ✓ Designed for full range flexibility and multi-positioning requirements
- ✓ Includes cushions, two 470-0300-0028 Clark Socket clamps and securing straps

Levitators Stirrups | 470-1330-0059

- ✓ Weight capacity is 600 lbs
- ✓ Includes cushions, clamps and straps
- ✓ Will never slip out of place once automatically locked into position
- ✓ Effortlessly position the patient quickly & efficiently with Gas Shock



Lithotomy Adult Leg Holders | 470-1330-0078

- ✓ Integral clamp allows for precise adjustments from 23" - 41"
- ✓ 8" arc allow for additional leg clearance
- ✓ Weight capacity is 60 lbs
- ✓ Complete with one pair of straps

Replacement Straps (pair) | 470-1330-0081



Remanufactured Allen Medical Yellofins® Stirrups with Allen Medical Lift-Assist® | 470-1330-0060-R

- ✓ Allows easy adjustment of abduction and lithotomy while maintaining the sterile field
- ✓ Includes cushions, straps & two 470-0300-0023 Side Rail/Foot Extension clamps

Pediatric Stirrup Set | 470-1330-0064

- ✓ The solution for pediatric laparoscopy and other pediatric procedures requiring lithotomy positioning, especially where intra-operative re-positioning may be required
- ✓ Complete with 5 interchangeable boot sizes, cushions and clamps
- ✓ PediStirrup represents the exclusive device designed for safe and flexible lithotomy positioning of babies (starting with newborn) and children with body lengths up to 1.40 meters



Maquet 1420/1425 Orthostar Leather Traction Boots | 470-1330-0001

Steris Orthovision/Orthographic Foot Traction Boots | 470-1330-0093



Yellofins® is a registered trademark of Allen Medical Systems, Inc.
Lift-Assist® is a registered trademark of Allen Medical Systems, Inc.

12 *Leg Supports*



Under Knee Crutch Supports (pair) | 470-1330-0077

- ✓ Premium soft pad reduces popliteal pressure
- ✓ Provides excellent patient comfort and precise surgical positioning for short duration procedures
- ✓ Complete with two 470-0300-0028 Clark Socket clamps

Arthroscopic Stress Post | 470-0930-0042

Arthroscopic Stress Post | 470-0930-0043

- ✓ Includes 470-0300-0028 Clark Socket clamp



Leg Prep Positioner

- ✓ Versatile leg prepper design, with 360° rotation
- ✓ Adjustable height of 12" - 27"
- ✓ Fits all surgical tables
- ✓ Requires 470-0300-0028 Clark Socket clamp

Standard Leg Prep Positioner with Pad | 470-1330-0082

- ✓ Leg cradle measures 8" H x 5" W x 3 ½" D

Bariatric Uni-Prepper Prep Unit with Pad | 470-1330-0083

- ✓ Leg cradle measures 8" H x 10.5" W x 3 ½" D



Universal Split Leg Accessory | 470-1330-0084

- ✓ Provides a safe system for positioning patients for laparoscopic, lower extremity orthopedic and vascular procedures
- ✓ Radiolucent, which allows both the surgeon and imaging equipment excellent access to the operative site

Steris I.A. Extender | 470-0400-0432

- ✓ Complete with 2" Sof-Care cushion



Steris 3080/3085 Patient Transfer Board | 470-0400-0433

- ✓ Complete with 2" Sof-Care cushion

Deluxe Arthroscopic Legholder System | 470-1330-0054

- ✓ Ratchet style strap easily adjusts for precise compression
- ✓ Holds legs up to 35" in circumference
- ✓ Requires one 470-0300-0028 Clark Socket clamp



14 *Positioning Devices*



Deluxe Lateral Brace Set | 470-0930-0044

- ✓ Fully adjustable vertically, horizontally and diagonally
- ✓ Includes three 470-0300-0028 Clark Socket clamps

Epidural Positioning Device | 470-0930-0057

- ✓ Includes positioning device and stand
- ✓ Allows positioning patients correctly and comfortably by encouraging cervical, thoracic and lumbar flexion
- ✓ Chest is supported, making the back immobile



Remanufactured OSI Wilson Spinal Frame | 470-0930-0059

Vacuum Activated Bean Bag Positioners

- ✓ Black polyurethane film is extremely durable
- ✓ Simply hook up to a suction line to evacuate air and product molds to shape
- ✓ Open valve and product returns to relaxed state
- ✓ Valve is replaceable if damaged

Rectangular design, size 20" x 20" | 470-0910-0002

Horseshoe design, size 30" x 40" | 470-0910-0003

Horseshoe design, size 46" x 46" | 470-0910-0005



Restraint Straps 15



Conductive Restraint Strap with Single Buckle | 470-1100-0026

Velcro Safety Restraint Strap (Pair) | 470-1100-0027



Conductive Restraint Strap with Velcro Lock | 470-1100-0028

Conductive Patient Restraint Strap | 470-1100-0030



16 *Restraint Straps*



Padded Restraint Strap with Dual Buckle | 470-1100-0032

Bariatric Padded Restraint Strap | 470-1100-0033



Armboard Restraint Straps (pair) | 470-1100-0034

Bariatric Nissen Restraint Strap Set | 470-1100-0036





Remanufactured Tenet T-Max Shoulder Chair | 470-0930-0047-R

- ✓ Designed to attach to the rails of almost any surgical table
- ✓ Superior head adjustment allows for infinite positions of tilt, rotation and forward slide of the head
- ✓ Sliding shoulder pad and posterior placement allows for unparalleled access to the surgical site
- ✓ Complete with two 470-0300-0021 side rail clamps

Remanufactured Steris Shoulder Chair Positioner | 470-0930-0058

- ✓ Removable shoulder panels
- ✓ Floating head support
- ✓ Attaches to any model of O.R. Table
- ✓ Includes cushions and clamps



18 *Storage Carts & Holders*



Lift Assist Stirrup Storage and Transport Cart | 470-1200-0018

- ✓ Six rails and a deep bin afford plenty of storing space
- ✓ Tip-resistant weighted base and reduced vibration during transport
- ✓ Polymer base is dent-resistant

Wall Mounted Stirrup Rack | 470-1200-0013

- ✓ Provides convenient storage for your stirrups, armboards, clamps, straps and other table accessories
- ✓ Dent-resistant polymer construction is easily cleanable
- ✓ Unique key slots accept $\frac{1}{4}$ " x 1" flat and $\frac{5}{8}$ " round mounting posts without the need for clamps





Flexible Urology Drain Hoop With Clamps | 470-0500-0020

- ✓ Provides excellent fluid containment for hysteroscopy or cysto procedures on any table
- ✓ Latex free drapes have a tissue retention filter and a flexible hose that is compatible with all suction systems for accurate fluid measurement
- ✓ Includes two 470-0300-0021 side rail clamps

Non-Sterile Replacement Disposable Drain Bags | 420-0000-0063
(case of 10)

Sterile Replacement Disposable Drain Bags | 420-0000-0064
(case of 20)



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.



Cat. No. 03-317