

LED EXAMINATION LIGHT WALL MOUNTED



Introduction

Skytron TR-7W examination light provides brilliant white light with just seven (7) high efficiency LEDs. The 12" (305mm) polymer lighthead provides an easy to position, exceptionally durable, fixture creating the perfect balance of performance, maneuverability, and dependability. Applications include Emergency, Labor and Delivery, General Examination, and Imaging Departments.

Standard Features

- Sealed lighthead provides smooth surface for easy wipedown between cases
- Lighthead mounted intensity control provides infinite range of intensity adjustments
- Spring balanced articulating arm offers easy positioning and stability
- Removable sterilizable handle
- Total reach of 68" (1727mm)


Accessories

Additional Sterilizable Handles (D51-040-000)

Architectural Specification

Lighting system should consist of a 12" (305mm) diameter LED exam light with sterilizable handle. Includes polymer lighthead with seven (7) high efficiency LEDs (55,000 Lux) and a minimum life of 30,000 hrs. Wall mounted, with a total reach of 68" (1727mm). 280° rotation of the head.

Technical Drawing



Technical Data

Electrical: 100-120 VAC or 200-240 VAC, 50/60 Hz, single phase. Power input connector is designed for 20 AWG (minimum) to 14 AWG (maximum) solid core wire only.

Intensity: 55,000 Lux at 1m

Light Field Diameter at 1m: 6-19/64" (160mm)

Color Temperature: 4500K

CRI: >95

LED Life: 30,000 hours

Total Reach: 78" (1980mm)

Fixture Weight: 46 lbs (20.9 kg)

Moment Load: 155 ft-lbs (210 Nm)



DISTRIBUTED BY
MEDITEK

Here for L.I.F.E.

1.800.567.8400 www.meditek.ca
1810 Ellice Avenue Winnipeg, MB R3H 0B7