

LED Examination Lights

The Skytron TR-7 Series brings revolutionary LED technology to exam and procedure lighting. With its advanced thermal management system, integrated

dimming controls, brilliant light output, and compact housing, the TR-7 meets the stringent requirements of both hospital designers and practitioners.


Glare-free and adjustable dimming from 100% to 10%

Removable sterilizable handle

Excellent color rendering CRI > 95

Ergonomic and aesthetic shape

Long lasting and energy efficient



TR-7W
Wall Mount



TR-7C
Ceiling Mount



TR-7F
Mobile Floor Stand

Features:

- LED lifetime 30,000 hours
- 55,000 lux / 1 m
- Fixed illuminated area 6.5" (160mm)
- No heat radiation
- Maintenance free
- Shadow and reflection free
- Available in ceiling or wall mount and mobile floor stand versions



Product Advantages

The TR-7 is ergonomically designed for comfort and ease of operation. Skytron's high quality components guarantee a durable, user-friendly light and LED technology allows for long periods of maintenance-free usage.

An ergonomically placed switch on the light head provides intensity and on/off control for the practitioner.



Ceiling mount with 70" (1575mm) reach allows light to be positioned for any exam.



Mobile floor stand offers excellent exam lighting that can be brought to any needed area with 38" (965mm) reach.

Wall mount with 68" (1727mm) reach for specific applications and low ceiling areas.



Homogenous illuminated area

Unique optic ensures optimum light field. The illuminated area meets the highest technical lighting requirements for exam lights. The reflection-, shadow- and glare-free lighting generates an ideal field of vision.



Color rendering

Accurate rendering of the colors of tissue - no color distortion.

- Ra > 95
- R9 > 90

This allows a clear recognition of different color planes within the tissue. The illuminated area appears natural and contrasted.

ORDERING INFORMATION	
Product Description	Model #
TR-7F, 120V 50/60HZ, Mobile Floor Stand	D15-174-200
TR-7W, 100V 50/60HZ Wall Mount	D15-175-200
TR-7C, 100-240V 50/60HZ, Ceiling Mount, 6" VST for 96"-102" Ceiling Height	D15-170-190
TR-7C, 100-240V 50/60HZ, Ceiling Mount, 12" VST for 103"-110" Ceiling Height	D15-170-200
TR-7C, 100-240V 50/60HZ, Ceiling Mount, 20" VST for 111"-118" Ceiling Height	D15-170-210
TR-7C, 100-240V 50/60HZ, Ceiling Mount, 28" VST for 119"-126" Ceiling Height	D15-170-220