

Transfer Boards



MEDITEK

Here for L.I.F.E.

1.800.567.8400 | info@meditek.ca | www.meditek.ca

Transfer Boards

Meditek transfer boards remove the need to lift when performing patient transfers. Our transfer boards are manufactured and designed so that staff do not have to use extra effort, which makes their job easier and safer.

Features and Benefits:

- ✓ High performance
- ✓ Radioluscent
- ✓ Non-toxic materials
- ✓ Anti-bacterially surface
- ✓ Safe and comfortable for patients
- ✓ Prevents back injuries to care workers
- ✓ Skin-friendly (no latex)
- ✓ Easy to clean and disinfect
- ✓ Disposable covers available
- ✓ Tear resistant nylon cloth
- ✓ Bridges gaps up to 20 cm
- ✓ Lightweight, easy suspendable



Transfer patients in 3 easy steps:

Turn the patient onto the side, slide the board under the body and return the patient to the original position.



Gently push the patient at the shoulder and hip, without applying force, and transfer.



Turn the patient onto the side. Remove the Rollbord starting at the feet.



Transfer Boards

Standard | 470-1400-0000

✓ Suitable for radiology, operating theatre, intensive care units and all special wards.

Size: 177 x 50 cm / 87 x 50 cm (70" x 20" / 35" x 20")

Weight: 2.7 kg



Bariatric | 470-1400-0001

✓ Suitable for Intensive Care (ICU) nursing services, on air and decubitus mattresses.

Size: 177 x 50 cm (70" x 20")

Weight: 3.8 kg



Mini-Gyn | 470-1400-0011

✓ Suitable for urology/gynaecology, wake-up stations and emergencies.

Size: 91 x 50 cm (36" x 20")

Weight: 1.4 kg



Surgiboard | 470-1400-0003

✓ Suitable for operating theatre, X-Ray, various wards (replaces aluminium rollers).

Size: 126 x 40 cm (50" x 16")

Weight: 1.5 kg



MRI/CT | 470-1400-0004

✓ Suitable for transfer of injured patients on/ from stretcher.

Size: 177 x 40 cm (70" x 16")

Weight: 2.5 kg



Leadership. Innovation. Focus. Excellence.



Cat. No. 08-915